

**E003:****Refrigerators and freezers**

PQS code:	E003/083
Type of appliance:	Icelined refrigerator
Manufacturer's reference:	GVR 225 AC
Manufactured in:	India
Company:	Godrej & Boyce MFG. Co. Ltd.
Address:	Plant 01, Pirojshanagar, Vikhroli West Mumbai 400 079 India
Telephone:	+91 22 67966201
Email:	medicalref@godrej.com ; hello@surechill.com
Web address:	http://www.godrejappliances./medical

Specifications

Climate zone:	Hot	Min rated ambient temp:	+10°C
Refrigerant:	R600a	Energy source:	Electric Mains 115/60Hz;#8;#Electric Mains 220- 240V 50/60Hz or Electric Mains 115/60Hz optional
Appliance tested at:	+43°C	Ext dimensions (HxLxD)	79.5 x 75.0 x 182.0 cm
Performance at:	+43°C	Fuel and cycle type:	Electric - absorption

REFRIGERATOR**FREEZER****Freeze protection: Grade A****Vaccine storage capacity:** 225 Litres
(Litres)**Gross volume (Litres):** 278 Litres**Holdover time (hours):** 55 hours**Energy consumption, stable running (kWh/24 hours):** 2.04 kWh/24 h**Energy consumption, cool down test (kWh/24 hours):** 2.38 kWh/24 h**Gross volume:****Waterpack freezing capacity:****Waterpack storage capacity:****Energy consumption during freezing:****Comments:**

"(i) Sure Chill is a patented technology licensed by The Sure Chill Company (UK)
(ii) Refrigerator operates on only 2.5 hrs/day power at 43°C ambient, post stabilization
(iii) Has a built in intelligent temperature control, early warning and energy management system
(iv) Unit requires addition of clean drinking water during installation, contact manufacturer for information"

Accessories:

"Included in price: 4 wire shelves ; 1 Padlock; 1 Funnel
Optional (charged extra): Extra wire shelves, plywood packaging, Voltage Stabilizer"

Spare parts (ref, price):

Compressor
Fan Motor
Electronic Controller Card
Fan Impeller
Voltage Stabilizer

Shipping volume:	1.3 m3	Shipping weight:	160 kg
Price year	2016	Incoterms	EXW
Price	Price 1 (1-10 Units) - USD 1800	Price 2 (10-99 Units) - USD 1765	Price 3 (100+ Units) - USD 1730
Quality standard:	- ISO 9001:2008 -	Specification ref.	
Verification report:	CE-JOB-DEL-16-000501- GVR 225AC	Verification laboratory:	Intertek India

Current PQS status: *pre-qualified: 21 Nov 2016***Validity until:** *May, 2017*

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.