



E003:

Refrigerators and freezers

PQS code:	E003/081
Type of appliance:	Icelined refrigerator
Manufacturer's reference:	GVR 75 Lite
Manufactured in:	India
Company:	Godrej & Boyce MFG. Co. Ltd.
Address:	Plant 01, Pirojshanagar, Vikhroli West Mumbai 400 079 India
Telephone:	+91 22 67966201
Email:	medicalref@godrej.com ; hello@surechill.com
Web address:	http://www.godrejappliances/medical

Specifications

Climate zone:	Hot	Min rated ambient temp:	+10°C
Refrigerant:	R600a	Energy source:	Electric Mains 115/60Hz;#1;#Electric Mains 220- 240V 50/60Hz
Appliance tested at:	+43°C	Ext dimensions (HxLxD)	151.5cm x 61.8cm x 77.4cm
Performance at:	+43°C	Fuel and cycle type:	Electric - compression

REFRIGERATOR

FREEZER

Freeze protection: Grade A

Vaccine storage capacity: (Litres)	72.5 litres
Gross volume (Litres):	99 litres
Holdover time (hours):	81 hrs
Energy consumption, stable running (kWh/24 hours):	1.47 kWh/24 h
Energy consumption, cool down test (kWh/24 hours):	2.77 kWh/24 h

Gross volume:	
Waterpack freezing capacity:	
Waterpack storage capacity:	
Energy consumption during freezing:	

- "(i) Sure Chill is a patented technology licensed by The Sure Chill Company (UK) guarantees never to freeze vaccines
- (ii) Data loggers records up to 30 days temperature reading with 5 customizable temperature alarm settings
- (iii) Has a built in intelligent temperature control, early warning and energy management system
- (iv) Unit requires addition of clean drinking water during installation,"

Comments:

Accessories:

"Included in price: 4 wire shelves ; 1 Padlock; 1 Funnel
Optional (charged extra): Extra wire shelves, plywood packaging, Voltage Stabilizer"

Spare parts (ref, price):

Included in the UNICEF Supply Division spare parts kit:
 1 No. Compressor - Part No. RCM 977 492 57
 1 No. Start Relay - Part No. RRE 977 492 58
 1 No. Condenser - Part No. RCN 977 571 06
 1 No. Electronic Control Card / thermostat - Part No. RCC 977 49260
 1 No. Voltage stabilizer - Part No. RCM 977 493 57
 1 No. Capacitor - Part No. RFN 977 492 63
 Other parts are available as spares on request.

Shipping volume:	0.86m3	Shipping weight:	115 Kg
Price year	2016	Incoterms	
Price	Price 1 (1-10 Units) - USD 972	Price 2 (10-99 Units) - USD 950	Price 3 (100+ Units) - USD 935
Quality standard:	- ISO 9001:2008 -	Specification ref.	
Verification report:	CE-JOB-DEL-16-000501-GVR 75 Lite	Verification laboratory:	Intertek India
Current PQS status:	<i>pre-qualified: 21 Nov 2016</i>	Validity until:	<i>May, 2022</i>

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.