



## E003:

## Refrigerators and freezers

<b>PQS code:</b>	<b>E003/082</b>
<b>Type of appliance:</b>	Icelined refrigerator
<b>Manufacturer's reference:</b>	GVR 99 Lite
<b>Manufactured in:</b>	India
<b>Company:</b>	Godrej & Boyce MFG. Co. Ltd.
<b>Address:</b>	Plant 01, Pirojshanagar, Vikhroli West Mumbai 400 079 India
<b>Telephone:</b>	+91 22 67966201
<b>Email:</b>	<a href="mailto:medicalref@godrej.com">medicalref@godrej.com</a> ; <a href="mailto:hello@surechill.com">hello@surechill.com</a>
<b>Web address:</b>	<a href="http://www.godrejappliances/medical">http://www.godrejappliances/medical</a>

## Specifications

<b>Climate zone:</b>	Hot	<b>Min rated ambient temp:</b>	+10°C
<b>Refrigerant:</b>	R600a	<b>Energy source:</b>	Electric Mains 115/60Hz;#1;#Electric Mains 220- 240V 50/60Hz
<b>Appliance tested at:</b>	+43°C	<b>Ext dimensions (HxLxD)</b>	170cm(H) x 61.8cm(W) x 77.4cm (D)
<b>Performance at:</b>	<b>+43°C</b>	<b>Fuel and cycle type:</b>	Electric - compression

### REFRIGERATOR

### FREEZER

**Freeze protection: Grade A****Vaccine storage capacity:** 98.5 L  
(Litres)**Gross volume (Litres):** 124 L**Holdover time (hours):** 59 Hrs 34 minutes**Energy consumption, stable running (kWh/24 hours):** 1.23 kWh/24 h**Energy consumption, cool down test (kWh/24 hours):** 1.73 kWh/24 h**Gross volume:****Waterpack freezing capacity:****Waterpack storage capacity:****Energy consumption during freezing:**

- "(i) Sure Chill is a patented technology licensed by The Sure Chill Company (UK) guarantees never to freeze vaccines
- (ii) Data loggers records up-to 30 days temperature reading with 5 customizable temperature alarm settings
- (iii) Has a built in intelligent temperature control, early warning and energy management system
- (iv) Unit requires addition of clean drinking water during installation,"

**Comments:**

**Accessories:**

"Included in price: 5 wire shelves ; 1 Padlock; 1 Funnel  
Optional (charged extra): Extra wire shelves, plywood packaging, Voltage Stabilizer"

**Spare parts (ref, price):**

Included in the UNICEF Supply Division spare parts kit:  
 1 No. Compressor - Part No. RCM 977 492 57  
 1 No. Start Relay - Part No. RRE 977 492 58  
 1 No. Condenser - Part No. RCN 977 571 07  
 1 No. Electronic Control Card / thermostat - Part No. RCC 977 49260  
 1 No. Voltage stabilizer - Part No. RCM 977 493 57  
 1 No. Capacitor - Part No. RFN 977 492 63

Other parts are available as spares on request. .

<b>Shipping volume:</b>	0.95m3	<b>Shipping weight:</b>	125 Kg
<b>Price year</b>	2016	<b>Incoterms</b>	
<b>Price</b>	Price 1 (1-10 Units) - USD 1050	Price 2 (10-99 Units) - USD 1030	Price 3 (100+ Units) - USD 1010
<b>Quality standard:</b>	- ISO 9001:2008 -	<b>Specification ref.</b>	
<b>Verification report:</b>	CE-JOB-DEL-16-000501-GVR 99 Lite	<b>Verification laboratory:</b>	Intertek India

**Current PQS status:**

*pre-qualified: 21 Nov 2016*

**Validity until:** *May, 2022*

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.